

Cloning	Estimated price
Simple cloning 1: Insert of one DNA fragment into a vector using sticky-end ligation. Insert should be provided within a cloning plasmid.	3500 SEK *
Simple cloning 2: Insert of a second DNA fragment into the above generated vector using sticky-end ligation. Insert should be provided within a plasmid.	3000 SEK *
Advanced cloning: Insert of DNA fragment where insert must be amplified using PCR primers from a cDNA library or plasmid. (includes Primer design and ordering).	4500 SEK *
Maxiprep of plasmid DNA (up to 1 mg)**	From 2600 SEK

General consultations associated with a regular order are free of charge. **Vector design consultations are 1000 SEK per hour.**

*Price for cloning includes sequence validation of the DNA fragment cloned into the vector and validation of backbone size via restriction enzyme digest. A quantity of miniprep product from final construct will be delivered.

**Prices include initial transformation of plasmid, verification of final products by restriction enzyme digestion.

Note: Please also be aware that occasionally cloning projects do not work as planned and alternative approaches must be used. In such special cases, we will be in contact with you to inform of the situation, including delays or changes in price.